



# Recombinant *Drosophila simulans* Protein maelstrom (mael)

<b>Product Code</b>	CSB-BP463559DMJ
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B4QL99
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Drosophila simulans</i> (Fruit fly)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAPKKQNGFM MFVNEWNRN AEGRRMTLAE AVYHCGTIWE KMDTQQRGPY NSDAKDANAA RRDKRGSLNG HGQVDKAQRE AAESLMDKAQ REAAESLMDM KRITTERLVLN AKMSHDLENA KFVFVAFNYF TKALTTDVYV PAEFAACEYS LKEGIRSIYS TMIDPGQIIF GQGSDALHNS STTHDLPLPP NALGEKNMVK LYRNILDYLT KCQGEGKTPI VFTPAENIGM VKSCFRYLEC EDDSRDGGGK IEVFDIQYLL FILKKEVMSV AGLNDEKINK FATDAFFKND FFEFTAGIAC QISSINTRDY HVATLIGLLP AQTINVYLG TLRSMHEVLS DNDTKLTGYI SFLFEVICGV ALMFWVLQKA RKELSETLLS ADYNNNEGKHP DVQV
<b>Source</b>	Baculovirus
<b>Target Names</b>	mael
<b>Protein Names</b>	Recommended name: Protein maelstrom
<b>Expression Region</b>	1-394
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.