



Recombinant Drosophila simulans Protein maelstrom (mael)

Product Code	CSB-MP463559DMJ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B4QL99
Product Type	Recombinant Protein
Immunogen Species	Drosophila simulans (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MAPKKQNGFM MFVNEWNRN AEGRRMTLAE AVYHCGTIWE KMDTQQRGPY NSDAKDANAA RRDKRGSLNG HGQVDKAQRE AAESLMDKAQ REAAESLMDM KRITTERLVLN AKMSHDLENA KFVFVAFNYF TKALTTDVYV PAEFAACEYS LKEGIRSIYS TMIDPGQIIF GQGSDALHNS STTHDLPLPP NALGEKNMVK LYRNILDYLT KCQGEGKTPI VFTPAENIGM VKSCFRYLEC EDDSRDGGGK IEVFDIQYLL FILKKEVMSV AGLNDEKINK FATDAFFKND FFEFTAGIAC QISSINTRDY HVATLIGLLP AQTINVYLGS TLRSMHEVLS DNDTKLTGYI SFLFEVICGV ALMFWVLQKA RKELSETLLS ADYNNNEGKHP DVQV
Source	Mammalian cell
Target Names	mael
Protein Names	Recommended name: Protein maelstrom
Expression Region	1-394
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.