



Recombinant *Saccharomyces cerevisiae* Enhancer of translation termination 1 (ETT1)

Product Code	CSB-MP465045SVP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B3LJA2
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain RM11-1a) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MAKRPLGLGK QSREKKRKVE SVEKKSDEPS RESTPVRSQM SVELDDADL DDELAQLKGL WSKYFHSDRD DEYVLNGIVH ECDRLRLSE EDKEIKKTLN DIFHGIYALA LSELTIFKAG DEEATEEKRK KDVSSFFENA IERVELGLSH FPESQFLKLV LAKIIFQRIP LEYISNLHLK SKDKKLDLVG QLEHGKKHFS IYENDTEFTF EILQMVNDLL DIVENFGREQ SIQEGIDSDN EEEEEIIDIE LEPEHPVYPL QQSLEANYEW LRNHFDKLLD NTNTDMKIYA SIANTLGELY LKKAEEPSKV FLSLQYDDGS SKKVSDKEAK NVQETALKHT KKALEYLEKA KLEDDPDTWV QVAEAYIDLG NLLDNESAEQ EEAYKTAEI LGKANKASHG KFQDVLDNFL QG
Source	Mammalian cell
Target Names	ETT1
Protein Names	Recommended name: Enhancer of translation termination 1
Expression Region	1-412
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.