



Recombinant *Drosophila ananassae* Probable cytosolic Fe-S cluster assembly factor GF22738 (GF22738)

| | |
|--------------------------|---|
| Product Code | CSB-EP465074DKY-B |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | B3N018 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Drosophila ananassae</i> (Fruit fly) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MSRFSGALQL TDIDDFITPA QECIKPIPIN KKNSKSGAKI SVQEDGYEE SGTGKQKLQK LEITLQDCLA CSGCITSAEG VLITQQSQEE LLKVLHENKK LKAEDNNLS RTVVFSVATQ PLISLAYRYK ISVEETARHL AGYFRSLGAD YVLSTKVADD LALLECRNEF LEKYRENNDL TMFSSSCPGW VCYAEKTHGN FILPHIATTR SPQQIMGVLV KQYLGKKLNL PASRIYHVTV MPCYDKKLEA SREDFESEVN ASRDVDCVIT SVEVEQMLLE VDQTLPAQEP ADLDWPWWDE RPELMVYSHE STYSGGYAEH VFRFAARELF NDARPIQLQI DQPRNRDLKE FFLEKDGIVL LRFAIANGFK NIQNLVQKLK RGKGSKFHFV EIMACPSGCI NGGAQVRPIS GQHVRELTQK LEELYTKLPR SQPDNIGTKK LYDDFFDFGFG TDKSHNLLHT SYHSVDKLIP ALNIKW |
| Source | E.coli |
| Target Names | GF22738 |
| Protein Names | Recommended name: Probable cytosolic Fe-S cluster assembly factor GF22738 |
| Expression Region | 1-476 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. |



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.