



Recombinant *Drosophila ananassae* Probable cytosolic Fe-S cluster assembly factor GF22738 (GF22738)

Product Code	CSB-MP465074DKY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B3N018
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila ananassae</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MSRFSGALQL TDIDDFITPA QECIKPIPIN KKNSKSGAKI SVQEDGYEYEE SGTGKQKLQK LEITLQDCLA CSGCITSAEG VLITQQSQEE LLKVLHENKK LKAEDNNLS RTVVFSVATQ PLISLAYRYK ISVEETARHL AGYFRSLGAD YVLSTKVADD LALLECRNEF LEKYRENNDL TMFSSSCPGW VCYAETHGN FILPHIATTR SPQQIMGVLV KQYLGKKLNL PASRIYHVTV MPCYDKKLEA SREDFSEVN ASRDVDCVIT SVEVEQMLLE VDQTLPAQEP ADLDWPWWDE RPELMVYSHE STYSGGYAEH VFRFAARELF NDARPIQLQI DQPRNRDLKE FFLEKDGIVL LRFAIANGFK NIQNLVQKLK RGKGSKFHFV EIMACPSGCI NGGAQVRPIS GQHVRELTQK LEELYTKLPR SQPDNIGTKK LYDDFFDFGFG TDKSHNLLHT SYHSVDKLIP ALNIKW
Source	Mammalian cell
Target Names	GF22738
Protein Names	Recommended name: Probable cytosolic Fe-S cluster assembly factor GF22738
Expression Region	1-476
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.