



Recombinant *Drosophila erecta* L-asparaginase-like protein GG20738 (GG20738)

Product Code	CSB-EP465079DLK-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B3NN96
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila erecta</i> (Fruit fly)
Purity	≥85% (SDS-PAGE)
Sequence	VWEKSHKY FRNRKLRERR IKWFGTKKTE IQPLLISTWN YTDANLQAWS VLQQGPQRTR MAVIQGCMAC QNQRCGRLLA GGSSPDSEGT LTLEAAIMDG ERLEYGAVAG MEGARNAILV ADAVLRHTKH SLLVGKSATK FARSLGYKEE FLTDGRTKNV WKKWRSNGCQ PNFWRDVRPS HTENCGPYTP LPEHLHQHPL HQEYAITQGQ HDQLAFLALD AEGKLHVASQ SSGAPFRIPG RVGDSAVPGA GIYADNKVGG AVASGDGDIL MRHLPAFLAV EAMRAGKKPD QAAEWVVQRL LRHNTEFNGA VVVVNRRGIY AAACAGLDEF HFVSSGGKEY LSMARVERIK CLEREDEVVD DGPKGLFPTI PEKQAVPRG
Source	E.coli
Target Names	GG20738
Protein Names	Recommended name: L-asparaginase-like protein GG20738
Expression Region	23-399
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.