



# Recombinant *Drosophila sechellia* Adenylosuccinate synthetase (GM23154)

<b>Product Code</b>	CSB-EP465442DMG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B4II68
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Drosophila sechellia</i> (Fruit fly)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSASATNGTH YEQLHQGRTK MYKSKVDVVL GAQWGDEGKG KVVDMLASDV DIVCRCQGGN NAGHTVVANG TEFDFHLLPS GVVNEKCVSV IGNGVVIHLP SLFDEVLKNE AKGLQHLENR LIISDRAHLV FDFHQHVDGM QEAEKGGKSL GTTKKGIGPA YSSKATRNGI RVGELLGDFN LFSEKFKSIV ATHVRLFPSI NVDVEAELAR YKDYADKVRP YVKDTICFLH TALRNGKTIL VEGANAAMLIDFGTYPYVT SSNCSIGGVL TGLGLPPQTI GEVIGVVKAY TTRVGDGPFPTTEQLNDIGDL LQTRGFEIGV TTKRKRRRCGW LDIPLLKYTS LVNGYTCICL TKLDILDITLP EIKVAVAYKK PNGEKLDHFP GTIAELGNIE VEYAVLPGWQ TSTEEVRNFK ELPENAQSYV RFLESELSVP VRWVGVGKGR ESIINVH
<b>Source</b>	<i>E.coli</i>
<b>Target Names</b>	GM23154
<b>Protein Names</b>	Recommended name: Adenylosuccinate synthetase Short name= AMPase Short name= AdSS EC= 6.3.4.4 Alternative name(s): IMP--aspartate ligase
<b>Expression Region</b>	1-447
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.