



Recombinant *Drosophila yakuba* Eukaryotic translation initiation factor 3 subunit G-2 (eIF3-S4-2)

Product Code	CSB-MP465498DMR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B4PP90
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila yakuba</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MKAFNTSWAD EVDADYVDGL PPSNEYIEGD YKYVTEYKFN DDGKKVKVVR TFKIEKQIVP KAVARRRSWV KFGDSRLDKP GPNSQTTMAS EEIFMQFIGS KEFDQTHETQ LDPGKNIACR RICNGEHWSV NCPYKGTSM SKTMMETKAN AAAAAISDP SKTGKYVPPF MKDGGGGPGG KNWGRGRERD DSSAVRISNL SESMTEADLE ELVKKIGPHT KMYLAREKNT GLCKGFAYVH FKFRQDAAAA IEVLNGHGYD HLILCWEWSK PQP
Source	Mammalian cell
Target Names	eIF3gb
Protein Names	Recommended name: Eukaryotic translation initiation factor 3 subunit G-2 Short name= eIF3g-2 Alternative name(s): Eukaryotic translation initiation factor 3 RNA-binding subunit 2 Short name= eIF-3 RNA-binding subunit-2 Eukaryotic
Expression Region	1-273
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.