



Recombinant Human Zinc finger protein 20 (ZNF20)

Product Code	CSB-BP026573HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P17024
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MMFQDSVAFE DVAVSFTQEE WALLDPSQKN LYRDVMQETF KNLTSVGKTW KVQNIIDEYK NPRRNLGLMR EKLCESEKESH HCGESFNQIA DDMLNRKTLT GITPCESSVC GEVGTGHSSL NTHIRADTGH KSSEYQEYGE NPYRNKECKK AFSYLDSFQS HDKACTKEKP YDGKECTETF ISHSCIQRHR VMHSGDGPYK CKFCGKAFYF LNLCLIHRI HTGVKPYKCK QCGKAFTRST TLPVHERHTH GVNADECKEC GNAFSFPSEI RRHKRSHTGE KPYECKQCGK VFISFSSIQY HKMHTHTGKPYE YECKQCGKAF RCGSHLQKHG RTHTGKPYE CRQCGKAFRC TSDLQRHEKT HTEDKPYGCK QCGKGFRCAS QLQIHERTHS GKPHCKEKC GKVFYFSSL RIHERHTHTGE KPHECKQCGK AFRYFSSLHI HERHTHTGDKP YECKVCGKAF TCSSSIRYHE RTHTGKPYE CKHCGKAFIS NYIRYHERTH TGEKPYQCKQ CGKAFIRASS CREHERHTHTI NR
Source	Baculovirus
Target Names	ZNF20
Protein Names	Recommended name: Zinc finger protein 20 Alternative name(s): Zinc finger protein KOX13
Expression Region	1-532
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.