



Recombinant Human Zinc finger protein 37A (ZNF37A)

Product Code	CSB-EP026704HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P17032
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MITSQGSVSF RDVTVGFTQE EWQHLDPAQR TLYRDVMLEN YSHLVSVGYC IPKPEVILKL EKGEEPWILE EKFPSQSHLE LINTSRNYSI MKFNEFNKGG KCFCDEKHEI IHSEEEPSEY NKNNGNSFWLN EDLIWHQKIK NWEQSFEYNE CGKAFPENSL FLVHKRGYTG QKTKYTEHG KTCDSFFIT HQTTHPRENH YGNECGENIF EESILLEHQS VYPFSQKLN L TPIQRTHSIN NIIYNECGT FFSEKLVHL QQRTHTGEKP YECHECGKTF TQKSAHTRHQ RTHTGGKPYE CHECGKTFYK NSDLIKHQRI HTGERPYGCH ECGKSFSEKS TLTQHQRTHT GEKPYECHEC GKTF SFKSVL TVHQKTHTGE KPYECYACGK AFLRKSDLIK HQRIHTGEKP YECNECGKSF SEKSTLTKHL RTHTGEKPYE CIQCGKFFCY YSGFTEHLRR HTGEKPFGCN ECGKTFRQKS ALIVHQRTHI RQKPYGCNQC GKSFCVKSKL IAHRHTHTGE KPYECNVCCK SFYVSKSLTV HQRIHLGRNP INVVNEGNYS G
Source	E.coli
Target Names	ZNF37A
Protein Names	Recommended name: Zinc finger protein 37A Alternative name(s): Zinc finger protein KOX21
Expression Region	1-561
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.