



# Recombinant Human Zinc finger protein 37A (ZNF37A)

<b>Product Code</b>	CSB-MP026704HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P17032
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MITSQGSVSF RDVTVGFTQE EWQHLDPAQR TLYRDVMLEN YSHLVSVGYC IPKPEVILKL EKGEWPWILE EKFPSQSHLE LINTSRNYSI MKFNEFNKGG KCFKDEKHEI IHSEEEPSEY NKNNGNSFWLN EDLIWHQKIK NWEQSFEYNE CGKAFPENSL FLVHKRGYTG QKTKYTEHG KTCDSFFIT HQTTHPRENH YGNECGENIF EESILLEHQS VYPFSQKLNLP TPIQRTHSIN NIIYNECGT FFSEKLVHL QQRTHTGEKP YECHECGKTF TQKSAHTRHQ RTHTGGKPYE CHECGKTFYK NSDLIKHQRI HTGERPYGCH ECGKSFSEKS TLTQHQRTHT GEKPYECHEC GKTFSEKSVL TVHQKTHTGE KPYECYACGK AFLRKSDLIK HQRIHTGEKP YECNECGKSF SEKSTLTKHL RTHTGEKPYE CIQCGKFFCY YSGFTEHLRR HTGEKPFGCN ECGKTFRQKS ALIVHQRTHI RQKPYGCNQC GKSFCVKSKL IAHHRHTHTGE KPYECNVCCK SFYVSKSLTV HQRIHLGRNP INVVNEGNYS G
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ZNF37A
<b>Protein Names</b>	Recommended name: Zinc finger protein 37A Alternative name(s): Zinc finger protein KOX21
<b>Expression Region</b>	1-561
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.