



# Recombinant Escherichia coli O157:H7 Glucose-1-phosphate adenylyltransferase (glgC)

<b>Product Code</b>	CSB-MP466794EOE
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B5YUI6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O157:H7 (strain EC4115 / EHEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MVSLEKNDHL MLARQLPLKS VALILAGGRG TRLKDLTNKR AKPAVHFGGK FRIIDFALSN CINSGIRRMG VITQYQSHTL VQHIQRGWSF FNEEMNEFVD LLPAQQRMKG ENWYRGTADA VTQNLDIIRR YKAEYVVILA GDHIYKQDYS RMLIDHVEKG ARCTVACMPV PIEEASAFGV MAVDENDKII EFVEKPANPP SMPNDPSKSL ASMGYVFDA DYLYELLEED DRDENSSSHDF GKDLIPKITE AGLAYAHPFP LSCVQSDPDA EPYWRDVGTL EAYWKANLDL ASVVPEDMY DRNWPIRTYN ESLPPAKFVQ DRSGSHGRTL NSLVSGGCVI SGSVVVQSVL FSRVRVNSFC NIDSAVLLPE VVVGRCRLR RCVIDRACVI PEGMVI GENA EEDARRFYRS EEGIVLV TRE MLRKLGHKQE R
<b>Source</b>	Mammalian cell
<b>Target Names</b>	glgC
<b>Protein Names</b>	Recommended name: Glucose-1-phosphate adenylyltransferase EC= 2.7.7.27 Alternative name(s): ADP-glucose pyrophosphorylase Short name= ADPGlc PPase ADP-glucose synthase
<b>Expression Region</b>	1-431
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.