



Recombinant Escherichia coli O157:H7 DNA ligase B (ligB)

Product Code	CSB-BP466801EOE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B5YWE6
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7 (strain EC4115 / EHEC)
Purity	>85% (SDS-PAGE)
Sequence	MKVWMAILIS ILCWQSSVWA VCPAWSPARA QEEIFRLQQQ IKQWDDDDYWK EGKSEVEDGV YDQLSARLTQ WQRCFVSEPR DVMMPLNGA VMHPVAHTGV RKMADKNALS LWMRERSDLW VQPKVDGVAV TLVYRDGKLN KAISRGNGLK GEDWTQKVSL ISAVPQTVSG PLANSTLQGE IFLQREGHIQ QQMGGINARA KVAGLMMRQD DSDTLNSLGV FVWAWPDGPQ LMTDRKELA TAGFTLTQRY TRAVKNADEV ARVRNEWKKA KLPFVTDGVV VRGAKEPESR HWLPGQAEWL VAWKYQPVAQ VAKVKAIQFA VGKSGKISVV ASLAPVMLDD KKVQRVNIGS VRRWQEWDIA PGDQILVSLA GQGIPRIDDV VWRGAERTKP TPPENRFNSL TCYFASDVCQ EQFISRLVWL GSKQVLGLDG IGEAGWRALH QTHRFEHIFS WLLLTPPEQLQ NTPGIAKSKS AQLWHRFNLA RKQPFTRWVM AMGIPLTRAA LNASDERSWS QLLFSTEQFW QQLPGTGSGR ARQVIEWKEN AQIKKLGSWL AAQQITGFEP
Source	Baculovirus
Target Names	ligB
Protein Names	Recommended name: DNA ligase B EC= 6.5.1.2 Alternative name(s): Polydeoxyribonucleotide synthase [NAD(+)] B
Expression Region	1-560
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.