



Recombinant Escherichia coli O157:H7 Anaerobic glycerol-3-phosphate dehydrogenase subunit B (glpB)

| | |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-EP466803EOE |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | B5YX30 |
| Product Type | Recombinant Protein |
| Immunogen Species | Escherichia coli O157:H7 (strain EC4115 / EHEC) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MRFDTIIMGG GLAGLLCGLQ LQKHGLRCAI VTRGQSALHF SSGSLDLLSH LPDGQAVTDI HSGLESLRQQ APAHPYSLG PQRVLDLACQ AQALIAESGA QLQGSVELPH QRITPLGTLR STWLSSPEVP VWPLPAKKIC VVGISGLMDF QAHLAAASLR ELDLAVETA E IELPELDVLR NNATEFRAVN IARFLDNEEN WPLLLDALIP VANTCEMILM PACFGLADDK LWRWLNEKLP CSLMLLPTLP PSVLGIRLQN QLQRQFVRQG GVWMPGDEVK KVTCKNGVVN EIWTRNHADI PLRPRFAVLA SGSFFSGGLV AERNGIREPI LGLDVLQTAT RGEWYKGDFF APQPWQQFGV TTDQTLRPSQ AGQTIENLFT IGSVLGGFDP IAQGCGGGVC AVSALHAAQQ IAQRAGGQQ |
| Source | E.coli |
| Target Names | glpB |
| Protein Names | Recommended name: Anaerobic glycerol-3-phosphate dehydrogenase subunit B Short name= Anaerobic G-3-P dehydrogenase subunit B Short name= Anaerobic G3Pdhase B EC= 1.1.5.3 |
| Expression Region | 1-419 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. |



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.