



Recombinant Escherichia coli O157:H7 Small heat shock protein ibpB (ibpB)

Product Code	CSB-MP466804EOE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B5YX92
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7 (strain EC4115 / EHEC)
Purity	>85% (SDS-PAGE)
Sequence	MRNFDLSPLM RQWIGFDKLA NALQNAGESQ SFPPYNIKES DDNHYRITLA LAGFRQEDLE IQLEGTRLSV KGTPEQPKEE KKWLHQGLMN QPFSLSFTLA ENMEVSGATF VNGLLHIDLI RNEPEPIAAQ RIAISERPAL NS
Source	Mammalian cell
Target Names	ibpB
Protein Names	Recommended name: Small heat shock protein ibpB Alternative name(s): 16 kDa heat shock protein B
Expression Region	1-142
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.