



Recombinant Escherichia coli O157:H7 Electron transport complex protein RnfB (rnfB)

Product Code	CSB-BP466850EOE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B5Z462
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7 (strain EC4115 / EHEC)
Purity	>85% (SDS-PAGE)
Sequence	MNAIWIAVAA VSLGLAFGA ILGYASRRFA VEDDPVVEKI DEILPQSQCG QCGYPGCRPY AEAISCNGEK INRCAPGGEA VMLKIAELLN VEPQLDGEA QELTPARMVA VIDENNCIGC TKCIQACPVD AIVGATRAMH TVMSDLCTGC NLCVDPCPTH CISLQPAET PDSWKWDLNT IPVRIIPVEH HA
Source	Baculovirus
Target Names	rsxB
Protein Names	Recommended name: Electron transport complex protein RnfB
Expression Region	1-192
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.