



Recombinant *Oryza sativa* subsp. *japonica* Putative beta-glucosidase 9 (BGLU9)

Product Code	CSB-MP467354OFG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7ECS8
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	GELPPISRR SFPKGFIFGT SSASYQCEGG AAEGGRGPSI WDTFTYQHPD KIADKSNQDV ADNTYHLYKE DVHMMKEMGM DAYRFSISWS RILPNGSLNG GVNIEGINYY NNLINELLLK GVQSFVTLFH YDTPQALEDK YNGFLSPNII NDYKDYAEIC FKEFGDRVKH WITFNEPWIF CSKAYASGTY APGRCSPWEM GKCSVGDSGR EPYTACHHQL LAHAETVRLY REKYQFTEEV VRQSQFIHDN DLHRRSAKLS FIIQNYLLLQ IHFQPGPGGR VCQYRH
Source	Mammalian cell
Target Names	BGLU9
Protein Names	Recommended name: Putative beta-glucosidase 9 Short name= Os4bglu9 EC= 3.2.1.21
Expression Region	32-316
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.