



# Recombinant Human CCAAT/enhancer-binding protein beta (CEBPB)

<b>Product Code</b>	CSB-MP005181HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P17676
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>MQRLVAWDPA CLPLPPPPPA FKSMEVANFY YEADCLAAAY GGKAAPAAPP AARPGPRPPA GELGSIGDHE RAIDFSPYLE PLGAPQAPAP ATATDTFEAA PPAPAPAPAS SGQHHDFLSD LFSDDYGGKN CKKPAEYGYV SLGRLGAAKG ALHPGCFAPL HPPPPPPPPP AELKAEPGFE PADCKRKEEA GAPGGGAGMA AGFPYALRAY LGYQAVPSGS SGSLSTSSSS SPPGTPSPAD AKAPPTACYA GAAPAPSQVK SKAKKTVDKH SDEYKIRRE RNNIAVRKSRD KAKMRNLETQ HKVLELTAEN ERLQKKVEQL SRELSTLRNL FKQLPEPLLA SSGHC</p>
<b>Source</b>	Mammalian cell
<b>Target Names</b>	CEBPB
<b>Protein Names</b>	Recommended name: CCAAT/enhancer-binding protein beta Short name= C/EBP beta Alternative name(s): Liver activator protein Nuclear factor NF-IL6 Transcription factor 5 Short name= TCF-5
<b>Expression Region</b>	1-345
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	The protein encoded by this intronless gene is a bZIP transcription factor which can bind as a homodimer to certain DNA regulatory regions. It can also form heterodimers with the related proteins CEBP-alpha, CEBP-delta, and CEBP-gamma. The encoded protein is important in the regulation of genes involved in immune and inflammatory responses and has been shown to bind to the IL-1 response element in the IL-6 gene, as well as to regulatory regions of several acute-phase and cytokine genes. In addition, the encoded protein can bind the promoter and upstream element and stimulate the expression of the collagen type I gene.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.