



# Recombinant Escherichia coli O157:H7 Multifunctional CCA protein (cca)

<b>Product Code</b>	CSB-MP469092EOE
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B5YR96
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O157:H7 (strain EC4115 / EHEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKIYLVGGAV RDALLGLPVK DRDWVVVGST PQEMLDAGYQ QVGRDFPVFL HPQTHEEYAL ARTERKSGSG YTGFTCYAAP DVTLEDDLKR RDLTINALAQ DDNGEIIDPY NGLGDLQNRL LRHVSPAFGE DPLRVLRVAR FAARYAHLGF RIADETLTLM REMTHAGELE HLTPERVWKE TENALTTRNP QVFFQVLRDC GALRVLFPEI DALFGVPAPA RWHPEIDTGI HTLMTLSMAA MLSPQVDVRF TTLCHDLGKG LTPPELWPRH HGHGPAGVKL VEQLCQLRV PNEIRDLARL VAEFHDLIHT FPMLNPKTIV KLFDSIDAWR KPQRVEQLAL TSEADVRGRT GFESADYPQG RWLREAWAVA QSVPTKAVVE AGFKGVEIRE ELTRRRIAAV ASWKEQRCPK PE
<b>Source</b>	Mammalian cell
<b>Target Names</b>	cca
<b>Protein Names</b>	Recommended name: Multifunctional CCA protein Including the following 4 domains: CCA-adding enzyme EC= 2.7.7.72 Alternative name(s): CCA tRNA nucleotidyltransferase tRNA CCA-pyrophosphorylase tRNA adenylyl-/cytidylyl-transfe
<b>Expression Region</b>	1-412
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.