



# Recombinant Escherichia coli O157:H7 Protein RecA (recA)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-MP469161EOE   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | B5Z2B0  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Escherichia coli O157:H7 (strain EC4115 / EHEC)   |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MAIDENKQKA LAAALGQIEK QFGKGSIMRL GEDRSMDVET ISTGSLSLDI<br>ALGAGGLPMG RIVEIYGPES SGKTTTLQV IAAAQREGKT CAFIDAEHAL<br>DPIYARKLGV DIDNLLCSQP DTGEQALEIC DALARSGAVD VIVVDSVAAL<br>TPKAEIEGEI GDSHMGLAAR MMSQAMRKLA GNLKQSNTLL IFINQIRMKI<br>GVMFGNPETT TGGNALKFYA SVRLDIRRIG AVKEGENVVG SETRVKVVKN<br>KIAAPFKQAE FQILYGEGIN FYGELVDLGV KEKLIKAGA WYSYKGEKIG<br>QGKANATAWL KDNPETAKEI EKKVRELLLS NPNSTPDFSV DDSEGVAETN<br>EDF |
| <b>Source</b>            | Mammalian cell  |
| <b>Target Names</b>      | recA  |
| <b>Protein Names</b>     | Recommended name: Protein RecA Alternative name(s): Recombinase A   |
| <b>Expression Region</b> | 1-353   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.             |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |