



# Recombinant Escherichia coli O157:H7 Probable cytosol aminopeptidase (pepA)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-YP469169EOE  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Uniprot No.</b>       | B5Z3L7   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | Escherichia coli O157:H7 (strain EC4115 / EHEC)  |
| <b>Purity</b>            | >85% (SDS-PAGE)  |
| <b>Sequence</b>          | MEFSVKSGSP EKQRSACIVV GVFEPRRLSP IAEQLDKISD GYISALLRRG<br>ELEGKPGQTL LLHHVNPVLS ERILLIGCGK ERELDERQYK QVIQKTINTL<br>NDTGSMEAVC FLTELHVKGR NNYWKVRQAV ETAKETLYSF DQLKTNKSEP<br>RRPLRKMVFN VPTRRELTSG ERAIQHGLAI AAGIKAAKDL GNMPPNICNA<br>AYLASQARQL ADSYSKNVIT RVIGEQQMKE LGMHSYLA VG QGSQNESLMS<br>VIEYKGNASE DARPIVLV GK GLTFDSGGIS IKPSEGMD EM KYDMCGAAAV<br>YGVMRMVAEL QLPINVIGVL AGCENMPGGR AYRPGDVLTT MSGQTVEVLN<br>TDAEGRLLVC DVLTYVERFE PEAVIDVATL TGACVIALGH HITGLMANHN<br>PLAHELIAAS EQSGDRAWRL PLGDEYQEQL ESNFADMANI GGRPGGAITA<br>GCFLSRFTRK YNWAHLDIAG TAWRSKGAKG ATGRPVALLA QFLLNRAGFN<br>GEE |
| <b>Source</b>            | Yeast  |
| <b>Target Names</b>      | pepA   |
| <b>Protein Names</b>     | Recommended name: Probable cytosol aminopeptidase EC= 3.4.11.1<br>Alternative name(s): Leucine aminopeptidase Short name= LAP EC= 3.4.11.10<br>Leucyl aminopeptidase   |
| <b>Expression Region</b> | 1-503  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | full length protein  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.  |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.   |



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.