



Recombinant *Borrelia burgdorferi* Flagellar assembly factor FliW (fliW)

Product Code	CSB-YP470048BUG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7J1B6
Product Type	Recombinant Protein
Immunogen Species	<i>Borrelia burgdorferi</i> (strain ZS7)
Purity	>85% (SDS-PAGE)
Sequence	MIDEKSIEFD FPEGILGFEN IKKFIKDSK YKPFSIMQSI NKDVSFLVTS PFNFLSEYLP NIQEKDWSDI KAKEEDEKVI LCIINMHVND YKDITANLKA PIIINKKKLL GKQAICTNEK YSLHHKVFKE
Source	Yeast
Target Names	fliW
Protein Names	Recommended name: Flagellar assembly factor FliW
Expression Region	1-130
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.