



Recombinant Escherichia coli O157:H7 23S rRNA (uracil (747)-C (5))-methyltransferase RlmC

Product Code	CSB-EP471453EOE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B5YSF1
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7 (strain EC4115 / EHEC)
Purity	>85% (SDS-PAGE)
Sequence	<pre>MQCALYDAGR CRSCQWITQP IPEQLSAKTA DLKNLLADFP VEEWCAPVSG PEQGFRNKAK MVSVGSVEKP LLGMLHRDGT PEDLCDCPLY PASFAPVFAA LKPFIARAGL TPYNVARKRG ELKYILLTES QSDGGMMLRF VLRSETKLAQ LRKALPWLQA QLPQLKVITV NIQPVHMAIM EGETEIYLTE QQALVERFND VPLWIRPQSF FQTNPAVASQ LYATARDWVR QLPVKHMWDL FCGVGGFGLH CATPDMQLTG IEIALEAIAC AKQSAELGL TRLQFQALDS TQFATAQGEV PELVLVNPPR RGIGKPLCDY LSTMAPRFII YSSCNAQTMA KDIRELPGYR IERVQLFDMF PHTAHYEVLTL LLVKM</pre>
Source	E.coli
Target Names	rlmC
Protein Names	Recommended name: 23S rRNA (uracil(747)-C(5))-methyltransferase RlmC EC= 2.1.1.189 Alternative name(s): 23S rRNA(m5U747)-methyltransferase
Expression Region	1-375
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.