



# Recombinant Zea mays Ninja-family protein 8

<b>Product Code</b>	CSB-BP471990ZAX
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B6TCE8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Zea mays (Maize)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MDDDNGLLELS LGLSLGGSSG KAKARDAPLE PKAEPQVEES SSKGGLQTPD APSGKFYQKS ADNHEQNSKQ RHSPVAPQFG NFWGQPGNSS GPVVDGSVEP VSHQPQLSSY QDGRMPINSG NNSEEQKPVS GNHKLPSSEEM NFQKKHQTAS DQPDAFSKSS DGGAKNAPIS ISTDDGSTGE NEDVAESEAE GSNSWLVAQR EDSAAGSVVN KGSDRKRSD AAADGFKGKR QPSFSGSESS SGKLTGPNLI SMQASNVVAL PYQVPGQVSG PPVLTNAPNF HPVCTVQLRP PTNGGLAVQT MGSASQVAFG YPAVQLPTLE TGSSWAFGAP PQALSSFTAK DKADQTGTKQ VDDGKKPQEA GASSAHAED EKKADRGLSL MGSAIRPGIA PNVKFGGSGS YPDLPWVSTI GAGPNGRTIS GVTYKFGRNE VKIVCACHGT HMSPEEFMRH ASADAPAQDN SATLPAFPAG NQATSAEN
<b>Source</b>	Baculovirus
<b>Protein Names</b>	Recommended name: Ninja-family protein 8
<b>Expression Region</b>	1-478
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.