



Recombinant *Saccharomyces cerevisiae* Altered inheritance of mitochondria protein 23, mitochondrial (AIM23)

Product Code	CSB-EP473583SVH
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B5VL46
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain AWRI1631) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	NASSTTDI FRNAMKRKRE LANLKEQSHG NVARNAAFPK EYIKRPKQVP RNATNRKKIL ITWSTGTDRA KEAANSVSE IFKKNHKGNI KVVDP TTHRI EPSNIRYFAK GIDLDKVGLS IVNVEQIDNE NQIPLVKIVE SRVALLKYSDFLAKKKEKEL MELGVLNKS YKNLVTDKED NLKHIKISWQ IESDDLKRQK AHEIVSLLKK GNKVTLYLDD KNNINSNNWL ENFEELDRSQ KGEPRLPES VFQKRAAVLE TLKEIVSEYA NDPVLLGNMN SKMIMKLI PK DVK PQNNDKR ALKELRKKER QEKLQKRIQR KKMNEM
Source	E.coli
Target Names	AIM23
Protein Names	Recommended name: Altered inheritance of mitochondria protein 23, mitochondrial
Expression Region	33-356
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.