



Recombinant Escherichia coli O157:H7 5-methyltetrahydropteroyltriglutamate--homocysteine methyltransferase (metE)

Product Code	CSB-MP473828EOE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B5YY78
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7 (strain EC4115 / EHEC)
Purity	>85% (SDS-PAGE)
Sequence	<p>MTILNHTLGF PRVGLRRELK KAQESYWAGN STRELLAVG RELRARHWDQ QKQAGIDLLP VGDFAWYDHV LTTSLLLGNV PARHQNKDGS VDIDTLFRIG RGRAPTGEPA AAAEMTKWFN TNYHYMVPEF VKGQQFKLTW TQLLDEVDEA LALGHKVKPV LLGPITWLWL GKVKGEQFDR LSLLNDILPV YQQVLAELAK RGIEWVQIDE PALVLELPQA WLDAYKPAYD ALQGQVKLLL TTYFEGVTPN LDTITALPVQ GLHVDLVHGK DDVVELHKRL PSDWLLSAGL INGRNVWRAD LTEKYVQIKD IVGKRDLWVA SSCSLLHSPI DLSVETRLDA EVKSWFAL QKCHELALLR DALNSGDtaa LAEWSAPIQA RRHSTRVHNP AVEKRLAAIT AQDSQRANVY EVRAEAQRAR FKLPAPWPTTT IGSFPQTTEI RTLRLDFKKG NLDANNYRTG IAEHIRQAIV EQERLGLDVL VHGEAERNM VEYFGEHLDG FVFTQNGWVQ SYGSRVCKPP IVIDVSRPA PITVEWAKYA QSMTDKPVKG MLTGPVTILC WSFPREDVSR ETIAKQIALA LRDEVADLEA AGIGIIQIDE PALREGLPLR RSDWDAYLQW GVEAFRINAA VAKDDTQIHT HMCYCEFNDI MDSIAALDAD VITIETSRSD MELLESEFEF DYPNEIGPGV YDIHSPNVPS VEWIEALLKK AAKRIPAERL WVNPDGGLKT RGWPETRAAL ANMVQAAQNL RRG</p>
Source	Mammalian cell
Target Names	metE
Protein Names	Recommended name: 5-methyltetrahydropteroyltriglutamate--homocysteine methyltransferase EC= 2.1.1.14 Alternative name(s): Cobalamin-independent methionine synthase Methionine synthase, vitamin-B12 independent isozyme
Expression Region	1-753
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.