



Recombinant Escherichia coli O157:H7 Pantothenate synthetase (panC)

Product Code	CSB-MP473838EOE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B5YZH0
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7 (strain EC4115 / EHEC)
Purity	>85% (SDS-PAGE)
Sequence	MLIETLPLL RQQIRRLRME GKRVALVPTM GNLHDGHEMVL VDEAKARADV VVVSIFVNPM QFDRPEDLAR YPRTLQEDCE KLNKRKVDLV FAPSVKEIYP NGTETHTYVD VPGLSTMLEG ASRPGHFRGV STIVSKLFNL VQPDIACFGE KDFQQLALIR KVMADMGFDI EIVGVPIMRA KDGLALSSRN GYLTAEQRKI APGLYKVLSS IADKLQAGER DLDEIIAAG QELNEKGFRA DDIQIRDADT LLEVSETSKR AVILVAAWLG DARLIDNKMV ELA
Source	Mammalian cell
Target Names	panC
Protein Names	Recommended name: Pantothenate synthetase Short name= PS EC= 6.3.2.1 Alternative name(s): Pantoate--beta-alanine ligase Pantoate-activating enzyme
Expression Region	1-283
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.