



# Recombinant Rat Mitochondrial-processing peptidase subunit alpha (Pmpca)

<b>Product Code</b>	CSB-EP018242RA-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P20069
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	SSGATYPN IPLSSPLPGV PKPIFATVDG QEKFETKVTT LDNGLRVASQ NKFGQFCTLG ILINSGSRYE AKYLSGIAHF LEKLAFSSTA RFDSKDEILL TLEKHGGICD CQTSRDITMY AVSADSKGLD TVVGLLADV LHPRLTDEEI EMTRMAVQFE LEDLNMRPDP EPLLTEMIHE AAFRENTVGL HRFCPVENIG KIDREVLHSY LKNYYTPDRM VLAGVGVEHE HLVECARKYL LGVQPAWGAP GAVWMLTAQW HSTRGGSSRW RETCQMSALR PPRFQSSHIY GGARELLLLLE EDFIPFAVLN MMMGGGGSFS AGGPGKGMFS RLYLNVLNRH HWMYNATSYH HSYEDTGLLC IHASADPRQV REMVEITKE FILMGRTVDL VELERAQTQL MSMLMMNLES RPFVIFEDVGR QVLATHSRKL PHELCTLRN VKPEDIKRNA SKMLRGKPAV AALGDLTDLP TYEHIQAALS SRDGRLPRTY RLFR
<b>Source</b>	E.coli
<b>Target Names</b>	Pmpca
<b>Protein Names</b>	Recommended name: Mitochondrial-processing peptidase subunit alpha EC=3.4.24.64 Alternative name(s): Alpha-MPP P-55
<b>Expression Region</b>	33-524
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.