



# Recombinant *Xenopus tropicalis* Wilms tumor protein homolog (wt1)

<b>Product Code</b>	CSB-MP474809XBF
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B5DE03
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Xenopus tropicalis</i> (Western clawed frog) ( <i>Silurana tropicalis</i> )
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGSDVRDMNA LLPPVSTLSG NSSCSMPVSS SGQWAPVLDL PPGAPYSSLP PHSFIKQEPT WNPDPHEDQC LSAFTVHFSG QFTGTAGACR YGAFGAPTPS QATTGQARMF PNSPYLSNCL DNQQGMRNQG YSAVAFDGTG SYGHTPSHHT AQFTNHSFKH EDPMGQQTSL GEQQYSVPPP VYGCHTPTDS CTGSQALLR TPYNSDNLYQ MTSQLECMW NQMNLGSSLK SHGTSYENDS HSTPMLYSCG GQYRIHTHGV FRGIQDVRRV PGVTPAIVRS STEANEKRPF MCAYPGCNKR YFKLSHLQMH SRKHTGEKPY QCDFKDCERR FSRSDQLKRH QRRHTGIKPF QCKTCQRKFS RSDHLKTHTR THTGEKPFSC RWPSCQKKFA RDELVRHHN MHQRNMTKLQ LAL
<b>Source</b>	Mammalian cell
<b>Target Names</b>	wt1
<b>Protein Names</b>	Recommended name: Wilms tumor protein homolog
<b>Expression Region</b>	1-413
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.