



Recombinant *Saccharomyces cerevisiae* Mitochondrial protein import protein ZIM17 (ZIM17)

Product Code	CSB-BP475332SVH
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B5VQB0
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain AWRI1631) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	NEV KKDDKKVHLG SFKVDKPKMM IAFTCKKCNT RSSHTMSKQA YEKGTVLISC PHCKVRHLIA DHLKIFHDHH VTVEQLMKAN GEQVSQDVGD LEFEDIPDSL KDVLGKYAKN NSENASQLPH PSQK
Source	Baculovirus
Target Names	ZIM17
Protein Names	Recommended name: Mitochondrial protein import protein ZIM17 Alternative name(s): Mitochondrial import inner membrane translocase subunit TIM15 mtHsp70 escort protein 1 mtHsp70-associated motor and chaperone protein TIM15/ZIM17 Short n
Expression Region	48-174
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.