



Recombinant Escherichia coli O157:H7 Chromosome initiation inhibitor (iciA)

Product Code	CSB-MP475530EOE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B5YQA9
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7 (strain EC4115 / EHEC)
Purity	>85% (SDS-PAGE)
Sequence	MKRPDYRTLQ ALDAVIRERG FERAAQKLCI TQSAVSQRIK QLENMFGQPL LVRTVPPRPT EQGQKLLALL RQVELLEEW LGDEQTGSTP LLLSLAVNAD SLATWLLPAL APVLADSPIR LNLQVEDETR TQERLRRGEV VGAVSIQHQA LPSCLVDKLG ALDYLFVSSK PFAEKYFPNG VTRSALLKAP VVAFDHLDDM HQAFLQQNFD LPPGSVPCHI VNSSEAFVQL ARQGTTCMI PHLQIEKELA SGELIDLTPG LFQRRMLYWH RFAPESRMMR KVTDALLDYG HKVLRQD
Source	Mammalian cell
Target Names	argP
Protein Names	Recommended name: Chromosome initiation inhibitor Alternative name(s): OriC replication inhibitor
Expression Region	1-297
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.