



Recombinant Escherichia coli O157:H7 D-tagatose-1,6-bisphosphate aldolase subunit kbaZ (kbaZ)

Product Code	CSB-YP475541EOE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B5YS22
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7 (strain EC4115 / EHEC)
Purity	>85% (SDS-PAGE)
Sequence	MKHLTEMVRQ HKAGKTNGIY AVCSAHPLVL EAAIRYASAN QTPLIEATS NQVDQFGGYT GMTPADFRGF VCQLADSLNF PQDALILGGD HLGPNRWQNL PAAQAMANAD DLIKSYVAAG FKKIHLDCSM SCQDDPIPLT DDIVAERAAR LAKVAEETCL EHFGEADLEY VIGTEVPVPG GAHETLSELA VTTPDAARAT LEAHRHAFEK QGLNAIWPRI IALVVQPGVE FDHTNVIDYQ PAKATALSQM VESYETLIFE AHSTDYQTPQ SLRQLVIDHF AILKVGPAIT FALREALFSL AAIEEELVPA KACSGLRQVL ENVMLNRPEY WQSHYHGDGN ARRLARGYSY SDRVRYYPWD SQIDDAFAHL VRNLADSPIP LPLISQYLPL QYVKVRSSEL QPTPRELIIN HIQDILAQYH TACEGQ
Source	Yeast
Target Names	kbaZ
Protein Names	Recommended name: D-tagatose-1,6-bisphosphate aldolase subunit kbaZ
Expression Region	1-426
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.