



Recombinant Escherichia coli O157:H7 UDP-4-amino-4-deoxy-L-arabinose--oxoglutarate aminotransferase (arnB)

Product Code	CSB-MP475570EOE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B5YX44
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7 (strain EC4115 / EHEC)
Purity	>85% (SDS-PAGE)
Sequence	MSEFLPFSRP AMGVEELAAV KEVLESGWIT TGPKNHAEQ AFCQLTGNQH AIAVSSATAG MHITLMALEI GKGDEVITPS LTWVSTLNMI SLQGATPVMV DVDRDTLMVT PEAIESAITP RTKAIIPVHY AGAPADINAI RAIGERYGIA VIEDAAHAVG TYYKGRHIGA KGTAIFSFA IKNITCAEGG LIVTDNENIA RQLRMLKFHG LGVDAYDRQT WGRAPQAEVL TPGYKYNLTD INAAIALTQL AKLEHLNTRR REIAQQYQQA LAALPFQPLS LPAWPHVHAW HLFIIIRVDEQ RCGISRDALM EALKERGIGT GLHFRAAHTQ KYRERFPSL SLPNTEWNSE RICSLPLFPD MTTADADRVI TALQQLAGQ
Source	Mammalian cell
Target Names	arnB
Protein Names	Recommended name: UDP-4-amino-4-deoxy-L-arabinose--oxoglutarate aminotransferase EC= 2.6.1.87 Alternative name(s): UDP-(beta-L-threo-pentapyranosyl-4"-ulose diphosphate) aminotransferase Short name= UDP-Ara4O aminotransferase UD
Expression Region	1-379
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.