



Recombinant Escherichia coli O157:H7 NADH-quinone oxidoreductase subunit B (nuoB)

Product Code	CSB-EP475577EOE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B5YXS7
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7 (strain EC4115 / EHEC)
Purity	>85% (SDS-PAGE)
Sequence	MDYTLTRIDP NGENDRYPLQ KQEIVTDPLE QEVNKNVFMG KLNDMVNWGR KNSIWPNYFG LSCCYVEMVT SFTAVHDVAR FGAEVL RASP RQADLMVVAG TCFTKMAPVI QRLYDQMLEP KVVISMGACA NSGGMYDIYS VVQGVDFKIP VDVYIPGCPP RPEAYMQALM LLQESIGKER RPLSWVVG DQ GVYRANMQSE RERKRGERIA VTNL RTPDEI
Source	E.coli
Target Names	nuoB
Protein Names	Recommended name: NADH-quinone oxidoreductase subunit B EC= 1.6.99.5 Alternative name(s): NADH dehydrogenase I subunit B NDH-1 subunit B
Expression Region	1-220
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.