



# Recombinant Escherichia coli Peptide methionine sulfoxide reductase MsrA (msrA)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP475831ENW   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | B6I2D2  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Escherichia coli (strain SE11)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MSLFDKKHLV SPADALPGRN TPMPVATLHA VNGHSMTNVP DGMEIAIFAM<br>GCFWGVVERLF WQLPGVYSTA AGYTGGYTPN PTYREVCSDG TGHAEAVRIV<br>YDPSVISYEQ LLQVFWENHD PAQGMRRQND QGTQYRSAIY PLTPEQDAAA<br>RASLERFQAA MLAADDDRRI TTEIANATPF YYAEDDHQQY LHKNPYGYCG<br>IGGIGVCLPP EA  |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | msrA  |
| <b>Protein Names</b>     | Recommended name: Peptide methionine sulfoxide reductase MsrA Short name= Protein-methionine-S-oxide reductase EC= 1.8.4.11 Alternative name(s): Peptide-methionine (S)-S-oxide reductase Short name= Peptide Met(O) reductase  |
| <b>Expression Region</b> | 1-212   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |