



# Recombinant Escherichia coli Anaerobic nitric oxide reductase transcription regulator norR (norR)

<b>Product Code</b>	CSB-EP475858ENW-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B6I695
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain SE11)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSFSVDVLAN IAIELQRGIG HQDRFQRLIT TLRQVLECD A SALLRYDSRQ FIPLAIDGLA KDVLGRRFAL EGHPRLEAIA RAGDVVRFPA DSELPDPYDG LIPGQESLKV HACVGLPLFA GQNLIGALTL DGMQPDQFDV FSDEELRLIA ALAAGALSNA LLIEQLESQN MLPGDAAPFE AVKQTQMIGL SPGMTQLKKE IEIVAASDLN VLISGETGTG KELVAKAIHE ASPRAVNPLI YLNCAALPES VAESELFHV KGAFTGAISN RSGKFEMADN GTLFLDEIGE LSLALQAKLL RVLQYGDIQR VGDDRSRVD VRVLAATNRD LREEVLGRF RADLFHRLSV FPLSVPLRE RGDDVILLAG YFCEQCRLRL GLSRVLSAG ARNLLQHYSF PGNVRELEHA IHRAVLSRA TRSGDEVILE AQHFAPPEVT LPPPEAAAVP VVKQNLREAT EAFQRETIRQ ALAQNHNNWA ACARMLETDV ANLHRLAKRL GLKD
<b>Source</b>	E.coli
<b>Target Names</b>	norR
<b>Protein Names</b>	Recommended name: Anaerobic nitric oxide reductase transcription regulator norR
<b>Expression Region</b>	1-504
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.