



Recombinant *Oryza sativa* subsp. *japonica* Peptide deformylase 1A, chloroplastic (PDF1A)

Product Code	CSB-BP476045OFG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B6RGY0
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	ASAGGGG WLSGLLGGKG GGGAPTAMTV TPGTVKAGDP VLHEPAQDVA PGDIPSEKVQ GVIDRMVAVM RKAPGVGLAA PQIGVPLKII VLEDTQEYIS YAPKKDIEAQ DRRPFDLLVI INPKLKTTSK RTALFFEGCL SVDGYRALVE RHLDVEVSGL DRNGRPIKVE ASGWQARILQ HECDHLEGTL YVDTMVPRTF RIVDNLDLPL PVGCPPIGAR
Source	Baculovirus
Target Names	PDF1A
Protein Names	Recommended name: Peptide deformylase 1A, chloroplastic Short name= OsPDF1A Short name= PDF 1A EC= 3.5.1.88
Expression Region	44-260
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.