



Recombinant *Saccharomyces cerevisiae* Increased recombination centers protein 19 (IRC19)

Product Code	CSB-MP477648SVH
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B5VMR4
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain AWRI1631) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MRKPSITITT AKAIITPDYT LIKSHSKYQL PSRFQKLDAD SPERSTVVKL FYRRFMLRKP FISNVKMKVD TYRDYVRYKF MKENYELKRY LVFNPDGLRS KIKVELLSNT KCCEKILPVT EMQRTLEFVL KSCSYLPETK AQKWDIARDN TYCRQILKNL LTMQYEKYRS ILHRGIGHDE LDVKFSLHKT TSSPLTKLNK TEKKKIPLFK VFSDFDTTLI YLNETLGTRL
Source	Mammalian cell
Target Names	IRC19
Protein Names	Recommended name: Increased recombination centers protein 19
Expression Region	1-230
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.