



Recombinant *Saccharomyces cerevisiae* rRNA biogenesis protein RRP36 (RRP36)

Product Code	CSB-YP477654SVH
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B5VSG4
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain AWRI1631) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MSYYFKNLKP DLNSDVEEDD GNLLSIMAN KSKREIDEQE SSDDELKTLS FGSLKKAETI IDEEDFKDTK PVHKKPITTT YREESFDEDD DSEDKSDEDA GFFEEDSEDE THHGQKVPKK KSKHAPIEQS SKKRVPVRN IPGLEIPRNK RSNLYQDIRF DKSTGKALDS SIIRKRYQFL DEYREKEIDE LQKLLQDRKF LSKIDQGERE EMEQRLKSMK SRLQSMKNKD LEREILKEYE SDMNKNNNTR YHLKKSEKRK VVQKWKFDHM KAKQREK VME RKRKKRLGKE FKQFEFHNRR
Source	Yeast
Target Names	RRP36
Protein Names	Recommended name: rRNA biogenesis protein RRP36 Alternative name(s): Ribosomal RNA-processing protein 36
Expression Region	1-300
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.