



# Recombinant Human Putative alpha-1-antitrypsin-related protein (SERPINA2)

<b>Product Code</b>	CSB-EP021058HU-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P20848
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	SLVEDPQGD AAQKTDTS HH DQGDWEDLAC QKISYNVTDL AFDLYKSWLI YHNQHVLVTP TSVAMAFRML SLGTKADTRT EILEGLNVNL TETPEAKIHE CFQQVLQALS RPDTRLQLTT GSSLFVNKSM KLVDTFLEDT KKLYHSEASS INFRDTEEAK EQINNYVEKR TGRKVVDLVK HLKKTSLAL VDYISFHGKW KDKFKAERIM VEGFHVDDKT IIRVPMINHL GRFDIHRDRE LSSWVLAQHY VGNATAFFIL PDPKMWQLE EKLTYSHLEN IQRAFDIRSI NLHFPKLSIS GTYKLRVPR NLGITKIFSN EADLSGVSQE APLKLSKAVH VAVLTIDEKG TEATGAPHLE EKAWSKYQTV MFNRPFLVII KEYITNFPLF IGKVVNPTQK
<b>Source</b>	E.coli
<b>Target Names</b>	SERPINA2
<b>Protein Names</b>	Recommended name: Putative alpha-1-antitrypsin-related protein Alternative name(s): Serpin A2
<b>Expression Region</b>	22-420
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.