



# Recombinant Escherichia coli O157:H7 Ribosomal RNA small subunit methyltransferase F (rsmF)

<b>Product Code</b>	CSB-MP477855EOE
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B5YQX9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O157:H7 (strain EC4115 / EHEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAQHTVYFPD AFLTQMREAM PSTLSFDDFL AACQRPLRRS IRVNTLKTSV ADFLRLTAPY GWTLTPIPWC EEGFWIERDD EDALPLGSTA EHLSGLFYIQ EASSMLPVAA LFADGNAPQR VMDVAAAPGS KTTQIAARMN NEGAILANEF SASRVKVLHA NISRCGISNV ALTHFDGRVF GAALPEMFDA ILLDAPCSGE GVVRKDPDAL KNWSPESNQE IAATQRELID SAFHALRSGG TLVYSTCTLN REENESVCLW LKETYPDAVE FLPLGDLFPG ANKALTEEGF LHVFPQIYDC EGFFVARLRK TQAIPALPTP KYKVGNFPPS PVKDREAAQI RQAAASVGLN WDGNLRLWQR DKEVWLFPPV IEALIGKVERF SRLGIKLAET HNKGYRWQHE AVIALASPDN ENAFELTPQE AEEWYRGRDV YPQAAPVADD VLVTFQHQPI GLAKRIGSRL KNSYPRELVR DGKLFTSNA
<b>Source</b>	Mammalian cell
<b>Target Names</b>	rsmF
<b>Protein Names</b>	Recommended name: Ribosomal RNA small subunit methyltransferase F EC=2.1.1.178 Alternative name(s): 16S rRNA m5C1407 methyltransferase rRNA (cytosine-C(5)-)-methyltransferase RsmF
<b>Expression Region</b>	1-479
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.