



# Recombinant Escherichia coli O157:H7 Ubiquinone/menaquinone biosynthesis methyltransferase ubiE (ubiE)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-MP477897EOE   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | B5YY82  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Escherichia coli O157:H7 (strain EC4115 / EHEC)   |
| <b>Purity</b>            | ≥85% (SDS-PAGE)   |
| <b>Sequence</b>          | MVDKSQETTH FGFQTVAKEQ KADMVAHV FH SVASKYDVMN DLMSFGIHRL<br>WKRFTIDCSG VRRGQTVLDL AGGTGDLTAK FSRLVGETGK VVLADINESM<br>LKMGREKLRN IGVIGNVEYV QANAEALPPF DNTFDCITIS FGLRNVTDKD<br>KALRSMYRVL KPGGRLLVLE FSKPIIEPLS KAYDAYSFHV LPRIGSLVAN<br>DADSYRYLAE SIRMHPDQDT LKAMMQDAGF ESVDYYNLTA GVALHRGYK<br>F   |
| <b>Source</b>            | Mammalian cell  |
| <b>Target Names</b>      | ubiE  |
| <b>Protein Names</b>     | Recommended name: Ubiquinone/menaquinone biosynthesis methyltransferase ubiE EC= 2.1.1.163 EC= 2.1.1.201 Alternative name(s): 2-methoxy-6-polyprenyl-1,4-benzoquinol methylase Demethylmenaquinone methyltransferase  |
| <b>Expression Region</b> | 1-251   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |