



Recombinant Escherichia coli O157:H7 GMP reductase (guaC)

Product Code	CSB-BP477905EOE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B5YZD9
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7 (strain EC4115 / EHEC)
Purity	≥85% (SDS-PAGE)
Sequence	MRIEDLKLKLG FKDVLRPKR STLKSRSDVE LERQFTFKHS GQSWSGVPPII AANMDTVGTF SMASALASFD ILTAVHKHYS VEEWQAFINN SSADVLKHVM VSTGTSADDF EKTKQILDNL PALNFVCIDV ANGYSEHFVQ FVAKAREAWP TKTICAGNVV TGEMCEELIL SGADIVKVGI GPGSVCTTRV KTGVGYPQLS AVIECADA AH GLGGMIISDG GCTTPGDVAK AFGGGADFVM LGGMLAGHEE SGGRIVEENG EKFMLFYGMS SESAMKRHVG GVAEYRAAEG KTVKLPLRGP VENTARDILG GLRSACTYVG ASRLKELTKR TTFIRVQEQE NRIFNNL
Source	Baculovirus
Target Names	guaC
Protein Names	Recommended name: GMP reductase EC= 1.7.1.7 Alternative name(s): Guanosine 5'-monophosphate oxidoreductase Short name= Guanosine monophosphate reductase
Expression Region	1-347
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.