



Recombinant Escherichia coli O157:H7 Acetyl-coenzyme A carboxylase carboxyl transferase subunit alpha (accA)

Product Code	CSB-EP477915EOE-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B5Z0G4
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7 (strain EC4115 / EHEC)
Purity	≥85% (SDS-PAGE)
Sequence	MSLNFLDFEQ PIAELEAKID SLTAVSRQDE KLDINIDEEV HRLREKSVEL TRKIFADLGA WQIAQLARHP QRPYTLDYVR LAFDEFDELA GDRAYADDKA IVGGIARLDG RPVMIIGHQK GREKKEKIRR NFGMPAPEGY RKALRLMQMA ERFKMPIITF IDTPGAYPGV GAEERGQSEA IARNLREMSR LGVPVVCTVI GEGGSGGALA IGVGDKVNML QYSTYSVISP EGCASILWKS ADKAPLAAEA MGIAPRLKE LKLIDSIPE PLGGAHRNPE AMAASLKAQL LADLADLDVL STEDLKNRRY QRLMSYGYA
Source	E.coli
Target Names	accA
Protein Names	Recommended name: Acetyl-coenzyme A carboxylase carboxyl transferase subunit alpha Short name= ACCase subunit alpha Short name= Acetyl-CoA carboxylase carboxyltransferase subunit alpha EC= 6.4.1.2
Expression Region	1-319
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.