



Recombinant Escherichia coli O157:H7 Chaperone protein hscA (hscA)

Product Code	CSB-EP477918EOE-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B5Z0Z9
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7 (strain EC4115 / EHEC)
Purity	>85% (SDS-PAGE)
Sequence	MALLQISEPG LSAAPHQRR L AAGIDLGTTN SLVATVRSQQ AETLADHEGR HLLPSVVHYQ QQGHSVGYDA RTNAALDTAN TISSVKRLMG RSLADIQQRY PHLPYQFQAS ENGLPMIETA AGLLNPVRVS ADILKALAAR ATEALAGELD GVVITVPAYF DDAQRQGTKD AARLAGLHVL RLLNEPTAAA IAYGLDSGQE GVIAVYDLGG GTFDISILRL SRGVFEVLAT GGDSALGGDD FDHLLADYIR EQAGIPDRSD NRVQRELLDA AIAAKIALSD ADSVTVNVAG WQGEISREQF NELIAPLVKR TLLACRRALK DAGVEADEVL EVVMVGGSTR VPLVRERVGE FFGRPPLTSI DPKVVAIGA AIQADILVGN KPDSEMLLLD VIPLSLGLET MGGLVEKVIP RNTTIPVARA QDFTTFKDGQ TAMSIHVMQG ERELVDQCRS LARFALRGIP ALPAGGAHIR VTFQVDADGL LSVTAMEKST GVEASIQVKP SYGLTDSEIA SMIKDSMSYA EQDVKARMLA EQKVEAARVL ESLHGALAAD AALLSAAERQ VIDDAAAHLS EVAQGDDVDA IEQAIKNVDK QTQDFAARRM DQSVRRALKG HSVDEV
Source	E.coli
Target Names	hscA
Protein Names	Recommended name: Chaperone protein hscA Alternative name(s): Hsc66
Expression Region	1-616
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.