



Recombinant Escherichia coli O157:H7 2-keto-4-pentenoate hydratase (mhpD)

Product Code	CSB-BP477931EOE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B5Z2Q5
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7 (strain EC4115 / EHEC)
Purity	≥85% (SDS-PAGE)
Sequence	MTKHTLEQLA ADLRRAAEQG EAIAPLRDLI GIDNAEAAYA IQHINVQYDV AQGRRVVGRK VGLTHPKVQQ QLGVDQPDFG TLFADMCYGD NEIIPFSRVL QPRIEAEIAL VLNRDLPATD ITFDELYNAI EWWLPAAEVV GSRIRDWSIQ FVDTVADNAS CGVYVIGGPA QRPAGLDLKN CAMKMTRNNE EVSSGRGSEC LGHPLNAAVW LARKMASLGE PLRAGDIILT GALGPMVAVN AGDRFEAHIE GIGSVAAAFS SAAPKGSLS
Source	Baculovirus
Target Names	mhpD
Protein Names	Recommended name: 2-keto-4-pentenoate hydratase EC= 4.2.1.80 Alternative name(s): 2-hydroxypentadienoic acid hydratase
Expression Region	1-269
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.