



Recombinant Escherichia coli Cytoplasmic trehalase (treF)

Product Code	CSB-BP478160ENW
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B6I385
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain SE11)
Purity	>85% (SDS-PAGE)
Sequence	MLNQKIQNPN PDELMIEVDL CYELDPYELK LDEMIEAEPE PEMIEGLPAS DALTPADRYL ELFEHVQSAK IFPDSKTFPD CAPKMDPLDI LIRYRKVRRH RFDLRLKFVE NHFWLPEVYS SEYVSDPQNS LKEHIDQLWP VLTREPQDHI PWSSLLALPQ SYIVPGGRFS ETYYWDSYFT MLGLAESGRE DLLKCMADNF AWMIENYGHI PNGNRYYLS RSQPPVFALM VELFEEDGVR GARRYLDHLK MEYAFWMDGA ESLIPNQAYR HVVRMPDGSL LNRYWDDRDT PRDESWLEDV ETAKHSGRPP NEVYRDLRAG AASGWDYSSR WLRDTGRLAS IRTTQFIPID LNAFLFKLES AIANISALKG EKETEALFRQ KASARRDAVN RYLWDDENGI YRDYDWRREQ LALFSAAIV PLYVGMANHE QADRLANAVR SLLTPGGIL ASEYETGEQW DKPNGWAPLQ WMAIQGFKMY GDDLLGDEIA RSWLKTVNQF YLEQHKMIEK YHIADGVPRE GGGGEYPLQD GFGWTNGVVR RLIGLYGEP
Source	Baculovirus
Target Names	treF
Protein Names	Recommended name: Cytoplasmic trehalase EC= 3.2.1.28 Alternative name(s): Alpha,alpha-trehalase Alpha,alpha-trehalose glucohydrolase
Expression Region	1-549
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.