



Recombinant Escherichia coli ADP-L-glycero-D-manno-heptose-6-epimerase (hldD)

Product Code	CSB-EP478162ENW
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B6I3J9
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain SE11)
Purity	>85% (SDS-PAGE)
Sequence	<p>MIIVTGGAGF IGSNIVKALN DKGITDILVV DNLKDGTKFV NLVDLNIADY MDKEDFLIQI MAGEEFGDVE AIFHEGACSS TTEWDGKYMM DNNYQYSKEL LHYCLEREIP FLYASSAATY GGRTSDFIES REYEKPLNVY GYSKFLFDEY VRQILPEANS QIVGFYFNV YGPREGHKGS MASVAFHLNT QLNNGESPKL FEGSENFKRD FVYVGDVADV NLWFLENGVS GIFNLGTGRA ESFQAVADAT LAYHKKGQIE YIPFPDKLKG RYQAFTQADL TNLRAAGYDK PFKTVAEGVT EYMAWLNLRDA</p>
Source	E.coli
Target Names	hldD
Protein Names	Recommended name: ADP-L-glycero-D-manno-heptose-6-epimerase EC=5.1.3.20 Alternative name(s): ADP-L-glycero-beta-D-manno-heptose-6-epimerase Short name= ADP-glyceromanno-heptose 6-epimerase Short name= ADP-hep 6-epimerase
Expression Region	1-310
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.