



Recombinant Chicken Adenylosuccinate lyase (ADSL)

| | |
|--------------------------|---|
| Product Code | CSB-MP001397CH |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P21265 |
| Product Type | Recombinant Protein |
| Immunogen Species | Gallus gallus (Chicken) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MATPCAEEDP LARYRSPLVS RYASAEMGFN FSERKKFGTW RRLWLYLAQA EKSLGLPITD ESIKEMEANL DNIDFKMAAE EEKLRHDVM AHVHTFAHCC PKAAAIHHLG ATSCYVGDNT DLIVLRDGFN LLLPKLARVI SRLADFAETH ADLPTLGFTH YQPAQLTTVG KRCCLWIQDL CMDLQNLERA RDDLRFRGVK GTTGTQASFL QLFEGDHSKV EELDRLVTAK AGFKRSYMTV GQTYSRKVDI EVLSVLASLG ASVHKICTDI RLLANLKEIE EPFEKDQIGS SAMPYKRNP RSERCCSLAR HLMTLVLDPL QTASVQWFER TLDDSANRRV CLAEAFITAD IILSTLQNIS EGLVVYPKVI DRRIRQELPF MATENIIMAM VKAGGNRQDC HEKIRVLSQQ AAVVKQEGG DNDFIARVRA DPYFSPHHEH LDSLLDPSSF TGRAPQQVAK FLKEEVRPAL IPYQSMGGK IELTL |
| Source | Mammalian cell |
| Target Names | ADSL |
| Protein Names | Recommended name: Adenylosuccinate lyase Short name= ASL EC= 4.3.2.2 Alternative name(s): Adenylosuccinase Short name= ASase |
| Expression Region | 1-485 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full length protein |
| Target Details | Adenylosuccinate lyase is involved in both de novo synthesis of purines and formation of adenosine monophosphate from inosine monophosphate. It catalyzes two reactions in AMP biosynthesis: the removal of a fumarate from succinylaminoimidazole carboxamide (SAICA) ribotide to give aminoimidazole carboxamide ribotide (AICA) and removal of fumarate from adenylosuccinate to give AMP. Adenylosuccinase deficiency results in succinylpurinemic autism, psychomotor retardation, and , in some cases, growth retardation associated with muscle wasting and epilepsy. Two transcript variants encoding different isoforms have been found for this gene. |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final |



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.