



# Recombinant Escherichia coli O157:H7 Glycine--tRNA ligase beta subunit (glyS)

<b>Product Code</b>	CSB-EP481184EOE-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B5YVL1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O157:H7 (strain EC4115 / EHEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSEKTFLVEI GTEELPPKAL RSLAESFAAN FTAELDNAGL AHGTVQWFAA PRRLALKVAN LAEAQPDREI EKRGPAIAQA FDAEGKPSKA AEGWARGCGI TVDQAERLTT DKGEWLLYRA HVKGESTEAL LPNMVATSLA KLPIPKLMRW GASDVHFVRP VHTVTLLLGD KVIPATILGI QSDRVIRGHR FMGEPEFTID NADQYPEILR ERGKVIADYE ERKAKIKADA EEAARKIGGN ADLSESLLEE VASLVEWPVV LTAKFEEKFL AVPSEALVYT MKGDQKYFPV YANDGKLLPN FIFVANIESK DPQQIISGNE KVVPRRLADA EFFFNTDRKK RLEDNLPRLQ TVLFQQQLGT LRDKTDRIQA LAGWIAEQIG ADVNHATRAG LLSKCDLMTN MVFEFTDTQG VMGMHYARHD GEAEVAVAL NEQYQPRFAG DDLPSNPVAC ALAIADKMDT LAGIFGIGQH PKGDKDPFAL RRAALGVLRI IVEKNLNDL QTLTEEAVRL YGDKLTNANV VDDVIDFMLG RFRAWYQDEG YTVDTIQAVL ARRPTRPADF DARMKAVSHF RTLEAAAALA AANKRVSNIL AKSDEVLSDR VNA STLKEPE EIKLAMQVVV LRDKLEPYFA EGRYQDALVE LAELREPVDA FFDKVMVMVD DKELRINRLT MLEKLRELFV RVADISLLQ
<b>Source</b>	E.coli
<b>Target Names</b>	glyS
<b>Protein Names</b>	Recommended name: Glycine--tRNA ligase beta subunit EC= 6.1.1.14 Alternative name(s): Glycyl-tRNA synthetase beta subunit Short name= GlyRS
<b>Expression Region</b>	1-689
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.